AdoptCell® -NK Kit (For research use only)

< MATERIALS PROVIDED >

COMPONENTS	VIALS	STORAGE CONDITIONS
AdoptCell® -NK-F	2	R.T.
AdoptCell® -NK-M	120 mL	4 °C
AdoptCell® -NK-S	1 vial	4 °C
AdoptCell® -NK-B	12 mL	R.T.

< MATERIALS & SUPPLIES REQUIRED >

Required:

- Human peripheral blood or Human peripheral blood mononuclear cells (PBMCs)
- Normal human AB-serum (KOHJIN BIO; 12181301)
- Ficoll-Paque™
- Phosphate buffered saline (PBS) (-)
- CD3⁺ cell-depletion kit (Dynabeads™ CD3 , Invitrogen; 11151D)
※ Please contact us if you use other products.
- Centrifuge tubes (15 mL, 50 mL conical tubes)
- Centrifuge
- CO₂ incubator
- Medical refrigerator (4 °C)
- Microscope
- Pipettes and pipette tips

Optional:

- OptiPrep™

< REAGENT PREPARATION >

1. Add 1 vial of AdoptCell®-NK-S to 120 mL of AdoptCell®-NK-M. (Store at 4°C)
2. Preparation before use : Add 5% normal human AB serum to the medium prepared in step 1 (henceforth referred to as final medium).

< METHODS >

▷ Day0

- ※ If you have PBMCs, start from step 2.
1. Isolate PBMCs from peripheral blood using Ficoll-Paque™ or OptiPrep™.
※ Please refer to OptiPrep™ Application Sheet C05.
 2. Deplete CD3⁺ cells from PBMCs. (According to the protocol of Dynabeads™ CD3, catalog no. 11151D)
 3. Adjust the density of CD3⁺ depleted PBMCs to 5×10^5 cells/mL using 12.0 to 15.6 mL of the final medium per an AdoptCell® -NK-F
 4. Incubate at 37 °C and 5% CO₂ humidified incubator.


▷ Day9

5. Add the final medium to each flask to a final volume of 50 mL. (If you start with 12 mL, add 38 mL of the final medium.)

▷ Day14

6. Harvest the cultured cells into the centrifuge tube by gently pipetting for several times.
7. Optional: When the cells are adhered, add 5 mL of AdoptCell®-NK-B per an AdoptCell® -NK-F and wait for 3 to 5 minutes at room temperature.
8. Optional: Tap the side of flask and pipette gently, and collect the cell suspension into the tube from step 6.
Note : Do not treat with AdoptCell®-NK-B for 5 minutes or more, and do not pipette vigorously.
9. Wash the AdoptCell® -NK-F with PBS and collect suspension into the tube from step 6.
10. Centrifuge at 500 g for 5 minutes at room temperature. Remove AdoptCell®-NK-B by replacing with a solution, such as PBS and use for any assays.


< Representative data >

Flow cytometric assay

PerCP/Cy5.5 anti-human CD3 antibody (BioLegend, 300430)
 PE/Cy7 anti-human CD16 antibody (BioLegend, 302016)
 Alexa Fluor®700 anti-human CD56 antibody (BioLegend, 318316)
 PerCP/Cy5.5 Mouse IgG1, κ isotype Ctrl antibody (BioLegend, 400150)
 PE/Cy7 Mouse IgG1, κ isotype Ctrl antibody (BioLegend, 400126)
 Alexa Fluor®700 Mouse IgG1, κ isotype Ctrl antibody (BioLegend, 400144)
 Analyzed using LSRFortessa (BD) with FlowJo Ver.10.4 software

Cytolytic Assay

Target cells: PKH26 (SIGMA, PHK26GL-1KT)

Death cells: 7-AAD

Analyzed using LSRFortessa (BD) with FlowJo Ver.10.4 software

%Lysis= [(% of target cell lysis - % of spontaneous cell death) / (% of maximum lysis - % of spontaneous cell death)] × 100
 spontaneous cell death : incubation of target cells in the cell alone
 maximum lysis : treatment of target cells with 10% formalin

< FREQUENTLY ASKED QUESTIONS >

- Q1. Why do cultured cells contain a large number of T cells?
 A1. Take care that the CD3⁺ cell proportion is less than 1% at the start of culture.
 If you use apheresis blood, we also recommend to remove CD34⁺ cells.
- Q2. Why do viability and activity appear to be low *in vivo*?
 A2. Although the stability of AdoptCell® NK has been confirmed in human blood, it has not been confirmed in mice, including NOG mice. Therefore, please try intraperitoneal administration.

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
PerCP/Cy5.5 anti-human CD3 Antibody (BioLegend, 300430)

PE/Cy7 anti-human CD16 Antibody (BioLegend, 302016)

Alexa Fluor®700 anti-human CD56 Antibody (BioLegend, 318316)


PerCP/Cy5.5 Mouse IgG1, κIsotype Ctrl Antibody (BioLegend, 400150)

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